HOW TO ORDER

ns sheet.	
meter that provides ease of mer, a tube tip cutter, and	
meter that provides ease of mer, a tube tip cutter, and	
owmeter that provides ease of mer, a tube tip cutter, and	
e specially calibrated flowmeter djustable wrench, a timer, a parately.	
Air flow indicator kit. Spot test system to determine the direction of air flow. Kit consists of aspirating bulb, carrying case and one box of No. 300 tubes.	
nd operation of the pump is	
in the container body.	
Kitagawa detector tube handbook.	
r op	

COMPRESSED BREATHING AIR TEST TUBES

600SP	Carbon monoxide - Breathing air, 5-100 ppm
601SP	Carbon dioxide - Breathing air, 100-3000 ppm
602SP	Oil mist - Breathing air, 0.3-5 mg/m ³
603SPA	Water vapor - Breathing air, 20-160 mg/m ³
604SP**	Oxygen - Breathing air, 2-24%

^{**}A 50 ml plastic syringe and a 1 meter vinyl tube are optionally necessary.

DETECTOR TUBES FOR DISSOLVED SUBSTANCES IN SOLUTION

10 tubes per box	
200SA	Sulphide ion, 2-1000 ppm
200SB	Sulphide ion, 0.5-10 ppm
201SA	Chloride ion, 10-2000 ppm
201SB	Chloride ion, 3-200 ppm
203S	Copper ion, 1-100 mg/L
204S	Cyanide ion, 0.2-5 ppm
205SL	Salinity, 0.01-0.8%
234SA	Free residual chlorine, 0.4-5 ppm



SGD-BK-CGA Compressed Breathing Air Sampling Kit